


Stroke	L1	L2
50	317.8	150
100	367.8	200
200	467.8	300
300	567.8	400

Notes: - Allowable static moment 22 Nm for all lengths
 - 'X' dimensions are user defined/variable



UNITS: mm [inch]

Unless otherwise specified, tolerance to be: CNS 4018 IT 14 Linear measures 0.5-3 ±0.1 3-6 ±0.1 6-30 ±0.2 30-120 ±0.3 120-315 ±0.5	NAME	DATE	DO NOT SCALE DRAWING	DIMENSIONS ARE METRIC
	DRAWN	M. Dank	May 25, 2012	TITLE: BSU-55- □ □ 0
	CHECKED			
PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MYOSTAT MOTION CONTROL INC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF COOLMUSCLE INC. IS PROHIBITED.	 Myostat Motion Control Inc. +1 905-836-4441 www.myostat.ca	SIZE B	PRT. NO. BSU-55- □ □ 0	REV X.0
SCALE: 1:2.2		WEIGHT: -	SHEET 1 OF 1	

